



Point-of-Care

User Guide

Showing Diagnosis View


SystemOne

Introduction


This guide will show you how DXS Point-of-Care reacts when you enter a diagnosis or the associated read code into your clinical system during a patient consultation.

Instructions


1. When entering a diagnosis or the associated Read code into your clinical system, the DXS Point-of-Care toolbar will expand and present links to relevant information in DXS Point-of-Care.




- In some cases, a banner might appear to highlight that important information is available.




- Click on the banner to access the information.
- The content within DXS Point-of-Care is filtered according to the diagnosis entered.



5. This is referred to as 'Diagnosis View'.



6. To view all of the content available in DXS Point-of-Care, click on the 'View All' button located at the bottom left corner of the screen.



Support



0800 028 0004 (option 1)



support-uk@dxs-systems.com



dxs-systems.co.uk

Training



0800 028 0004 (option 2)



training@dxs-systems.co.uk



training.dxs-systems.com

