



Point-of-Care

User Guide


How to Search in the Professional
Section

Introduction


This guide will show you how to search in the 'Professional' section in DXS Point-of-Care using the search bar.

Instructions


1. Open DXS Point-of-Care and navigate to the 'Professional' section.




2. In the 'Professional' section, the search bar allows you to search for articles based on their title or keywords.




3. Click on the search bar and type in a title or keyword.
4. The title or keyword can consist of single, multiple or partial words.




5. Adding an asterisk (*) before or after a word will search for any documents containing that word.



6. Press enter on your keyboard or click on the 'Find' button to initiate the search.




7. A list of search results will be displayed.




8. Results are displayed in categories denoting the type of document.


9. The categories, as well as the listed search results under each category, are displayed in alphanumerical order.




10. Any Favourite or Recently Accessed documents will be prioritised to the top of the list, making it quicker and easier to access them.




11. To filter your search results, select the chevron icon.



12. Filter buttons can be activated or deactivated by selecting or deselecting the tick box on the left of the filter button.
13. Please note that any changes to the filter selection or deselection will be retained and applied to the next search.



14. If you would like to refine your search, you can search for example, **'cancer +NICE'**.
15. The '+' sign indicates that the next word must be in the name of the publisher of the article.



16. This will only return articles related to cancer that are published by NICE.

The screenshot shows a professional software interface with a search bar at the top containing the text "cancer +NICE". Below the search bar, the results are displayed under the heading "Guidance". The results list various NICE guidance documents related to cancer treatment, such as "Abemaciclib With An Aromatase Inhibitor For Previously Untreated, Hormone Receptor-Positive, Her2-Negative, Locally...", "Abiraterone Treat Metastatic Hormone-Relapsed Prostate Cancer Before Chemo(TA387)", and "Advanced Breast Cancer - Diagnosis And Treatment (CG81)". The interface includes a navigation menu on the left with categories like "My Favourites", "Recent", "Practice Content", "DXS Media", "Accident and Emergency", "Alerts", "Allergy and Immunology", "Cancer", "Cardiovascular", "Cardiovascular Disease", "Carers", "Community Services", "Coronavirus", "Dental and Oral Health", "Dermatology", "Diabetes", "Diagnostics", "Diet and Nutrition", "Elderly Care", "End of Life and Palliative Care", "Endocrine disorders", "Endocrinology", "ENT", "Fertility", "Forms and Protocols", "Gastroenterology", "Gastrointestinal Disorders", "Gender Identity", "General", "General Health and Wellbeing", "Genetic Disorders", and "Genitourinary Disorders". The bottom of the screen shows a Windows taskbar with various application icons and a system tray displaying the time as 11:51 on 22/06/2020.

Support



0800 028 0004 (option 1)



support-uk@dxs-systems.com



dxs-systems.co.uk

Training



0800 028 0004 (option 2)



training@dxs-systems.co.uk



training.dxs-systems.com

