



Point-of-Care

Clinical Decision Support

DXS PoC Release Notes

Version 5.4.8 – Limited External
Release

Date: 6 October 2021

Document Owner: Zayne Glover

Content

| | |
|--|---|
| 1. Product Description | 2 |
| 2. Releases | 2 |
| 3. Notes on Deployment | 2 |
| 4. System Requirements | 3 |
| 4.1 Recommended Specification | 3 |
| 4.2 Remote Database Server Specification | 3 |
| 5. PoC 5.4.8 – Limited External Release | 4 |
| 5.1 Introduction | 5 |
| 5.2 Enhancements..... | 5 |
| 5.3 Fixes | 5 |
| 5.4 Document changes | 5 |
| 6. PoC 5.4.7 – Previous General Release | 6 |
| 6.1 Introduction | 7 |
| 6.2 Enhancements..... | 7 |
| 6.3 Fixes | 7 |
| 6.4 Known issues..... | 7 |
| 6.5 Document changes | 8 |

1. Product Description

The DXS Point-of-Care™ product is a content aggregator, which integrates with Electronic Patient Record (EPR) software, ensuring a prime position on the healthcare professional's desktop. This allows content providers to have a direct channel to the healthcare provider during the patient consultations. The DXS product enriches the content with its own relationship-linking process, part of which incorporates the appropriate country-specific diagnostic coding system.

The diagnosis is entered into the EPR and information in DXS Point-of-Care™ relating to that diagnosis is automatically populated for the healthcare provider. The healthcare provider is presented with the specific information related to the patient he or she is consulting, without having to search different databases, which expedites the review.

The DXS Point-of-Care™ product contains clinical data on diseases and procedures as well as patient information on drugs and diseases.

2. Releases

| Status | Version | Release Date | Database Schema |
|------------------|----------|----------------|-----------------------|
| Current Release | 5.4.8.36 | 8 October 2021 | Version 8 Revision 14 |
| Previous Release | 5.4.7.25 | 31 August 2021 | Version 8 Revision 14 |

3. Notes on Deployment

The upgraded files will be released electronically and will be downloaded with the daily content update. After the update, the upgrade will be pending and the application will need to restart to deploy the new files.

4. System Requirements

4.1 Recommended Specification

| DXS Point-of-Care | Recommended Requirements |
|--------------------------------|--|
| Processor | Intel i3 or better |
| RAM | 2GB RAM |
| Workstation – HDD – Free Space | 200MB |
| Server – HDD – Free Space | 4GB (Growth detailed in Operating Environment Guidance document, paragraph 2.5) |
| Resolution | 1280 × 720 |
| Display | 32-bit Colour |
| Operating System | Windows 7, 8, 8.1, 10 supported |
| Internet Connectivity | Required (HTTP access to uk.dxs-systems.com. HTTP GET and POST Methods allowed) |
| Internet Speed | 5MB (Download), 2MB (Upload) – Please see below for practices that have branch sites |

4.2 Remote Database Server Specification

The below Internet speeds are our minimum requirements if a branch site connects remotely via a WAN to the same database as the main site.

| DXS Point-of-Care | Recommended Requirements |
|--------------------------------|---|
| Processor | Intel i3 or better |
| RAM | 2GB RAM |
| Workstation – HDD – Free Space | 200MB |
| Server – HDD – Free Space | 4GB (Growth detailed in Operating Environment Guidance document, paragraph 2.5) |
| Resolution | 1280 × 720 |
| Display | 16 -bit Colour |
| Operating System | Windows 7, 8, 8.1, 10 supported |
| Internet Connectivity | Required (HTTP access to uk.dxs-systems.com. HTTP GET and POST Methods allowed) |
| Internet Speed | 10MB (Download), 10MB (Upload) |

5. PoC 5.4.8 – Limited External Release

| Application | Version | Date | Status |
|----------------|----------|------------|---------|
| DXS.exe | 5.2.8.36 | 01.10.2021 | Updated |
| DXSSch.exe | 5.0.1.22 | 29.10.2016 | – |
| DXScfg.exe | 2.0.0.25 | 25.02.2005 | – |
| DXSBar.exe | N/A | 14.06.2005 | – |
| Launch.exe | 2.2.7.25 | 01.10.2021 | Updated |
| APIProxy.exe | 2.4.7.25 | 01.10.2021 | Updated |
| RegLib.exe | N/A | 01.10.2021 | Updated |
| DXS Loader.exe | 1.0.0.44 | 17.03.2021 | – |
| Linker.exe | 1.1.6.93 | 01.10.2021 | Updated |
| HWUpSvc | 1.1.0.38 | 10.02.2021 | – |

| Components | Version | Date | Status |
|-------------------------------------|--------------|------------|------------------|
| Hwise32x.dll | 5.4.8.36 | 01.10.2021 | Updated |
| Hwise32x.eng | 5.4.8.36 | 01.10.2021 | Updated |
| DXSBar32.dll | N/A | | – |
| DCINPS.dxc | 3.2.8.36 | 01.10.2021 | Updated (Vision) |
| DCEMISWb.dxc | 2.2.8.36 | 01.10.2021 | Updated (EMIS) |
| DCTPPSO.dxc | 2.2.8.36 | 01.10.2021 | Updated (TPP) |
| Stats.dxp | 2.0.2.11 | 13.08.2013 | – |
| PDPL.dxp | 3.0.1.19 | 11.12.2018 | – |
| Cardiac.dxp | 1.0.2.16 | 10.12.2008 | – |
| ssleay32.dll | 1.0.2.15 | 11.12.2018 | – |
| libeay32.dll | 1.0.2.15 | 11.12.2018 | – |
| FirebirdSql.Data.FirebirdClient.dll | 4.6.1.0 | 06.01.2020 | – |
| HWInput.eng | 2.2.8.36 | 01.10.2021 | Updated |
| LatticeLog.dll | 1.4.1.1 | 06.01.2020 | – |
| LatticeLogBI.dll | 1.0.3.1 | 06.01.2020 | – |
| LatticeLogDAL.dll | 1.0.3.1 | 06.01.2020 | – |
| Microsoft.ServiceBus.dll | 3.0.51272.20 | 06.01.2020 | – |
| Newtonsoft.Json.dll | 10.0.1.20720 | 06.01.2020 | – |
| Lwxspell.dll | 1.3.1.0 | 27.02.2020 | – |
| Editorial.dxp | 1.0.2.67 | 29.07.2021 | – |

| Help Files | Date |
|------------|------------|
| DXS.chm | 17.05.2007 |
| PDPL.chm | 16.07.2009 |

5.1 Introduction

Version 5.4.8 provides a significant performance boost around saving forms to the patient record and introduces placeholder fields, allowing dynamic switching between the type of field that appears on a form based on data in the patient record.

5.2 Enhancements

1. Saving forms to the patient record is now completely done in the background.
2. Placeholder fields have been introduced.
3. Error detection and logging have been extended to include errors around saving to the patient record.
4. Ability added for DXS Point-of-Care to force users to use SNOMED at a predetermined date.

5.3 Fixes

1. TFS-4454 – Fix for the Email button not opening any mail client.

5.4 Document changes

| Category | Document | Updated |
|----------------------------|---|---------|
| Deployment | DXS Centrally Hosted SQL Server | |
| | DXS Operating Environment Guidance | |
| | DXS PoC Deployment Detail | |
| | DXS Practice Content Migration | |
| | Technical Survey Overview | |
| | TIO Questionnaire | |
| Knowledge Base | DXS PoC Migration – CSU Instructions | |
| | DXS Program Upgrades and Environment Overview | |
| Technical Overview | PoC Feedback Messaging Process | |
| System Requirements | Release Notes | X |
| | DXS Centrally Hosted SQL Server | |
| | DXS Operating Environment Guidance | |
| | DXS PoC Deployment Detail | |
| | Technical Survey Overview | |
| | TIO Questionnaire | |

6. PoC 5.4.7 – Previous General Release

| Application | Version | Date | Status |
|----------------|----------|------------|---------|
| DXS.exe | 5.2.7.25 | 31.08.2021 | Updated |
| DXSSch.exe | 5.0.1.22 | 29.10.2016 | – |
| DXScfg.exe | 2.0.0.25 | 25.02.2005 | – |
| DXSBar.exe | N/A | 14.06.2005 | – |
| Launch.exe | 2.2.7.25 | 31.08.2021 | Updated |
| APIProxy.exe | 2.4.7.25 | 31.08.2021 | Updated |
| RegLib.exe | N/A | 31.08.2021 | Updated |
| DXS Loader.exe | 1.0.0.44 | 17.03.2021 | – |
| Linker.exe | 1.1.6.93 | 31.08.2021 | Updated |
| HWUpSvc | 1.1.0.38 | 10.02.2021 | – |

| Components | Version | Date | Status |
|-------------------------------------|--------------|------------|------------------|
| Hwise32x.dll | 5.4.7.25 | 31.08.2021 | Updated |
| Hwise32x.eng | 5.4.7.25 | 31.08.2021 | Updated |
| DXSBar32.dll | N/A | | – |
| DCINPS.dxc | 3.2.7.25 | 31.08.2021 | Updated (Vision) |
| DCEMISWb.dxc | 2.2.7.25 | 31.08.2021 | Updated (EMIS) |
| DCTPSO.dxc | 2.2.7.25 | 31.08.2021 | Updated (TPP) |
| Stats.dxp | 2.0.2.11 | 13.08.2013 | – |
| PDPL.dxp | 3.0.1.19 | 11.12.2018 | – |
| Cardiac.dxp | 1.0.2.16 | 10.12.2008 | – |
| ssleay32.dll | 1.0.2.15 | 11.12.2018 | – |
| libeay32.dll | 1.0.2.15 | 11.12.2018 | – |
| FirebirdSql.Data.FirebirdClient.dll | 4.6.1.0 | 06.01.2020 | – |
| HWInput.eng | 2.2.7.25 | 31.08.2021 | Updated |
| LatticeLog.dll | 1.4.1.1 | 06.01.2020 | – |
| LatticeLogBI.dll | 1.0.3.1 | 06.01.2020 | – |
| LatticeLogDAL.dll | 1.0.3.1 | 06.01.2020 | – |
| Microsoft.ServiceBus.dll | 3.0.51272.20 | 06.01.2020 | – |
| Newtonsoft.Json.dll | 10.0.1.20720 | 06.01.2020 | – |
| Lwxspell.dll | 1.3.1.0 | 27.02.2020 | – |
| Editorial.dxp | 1.0.2.67 | 29.07.2021 | – |

| Help Files | Date |
|------------|------------|
| DXS.chm | 17.05.2007 |
| PDPL.chm | 16.07.2009 |

6.1 Introduction

Version 5.4.7 adds more error detection and logging, covering a wider product area, including email and updating errors. This version also includes various bug fixes.

6.2 Enhancements

1. Error detection and logging have been extended to include email, updating, and user-not-registered errors
2. The CCG approval date has been added to content.

6.3 Fixes

1. TFS-4421 – Fix for content of pathway blocks disappearing after an update
2. TFS-4431 – Fix for the Document History tab giving an access violation error.
3. TFS-4432 – Fix for BMI search fields not populating for Vision.
4. TFS-4434 – Fix for input errors not being cleared if the value is changed via another field.
5. TFS-4455 – Fix for the medications being sent to ExpertCare without DM+D codes when using Vision.

6.4 Known issues

1. PRB-18600 – PDP Appraisal Summary values slightly difficult to read when using the 'High Contrast' (non-default) Windows theme.
2. PRB-18711 – Content being cut off from the right side in print. Only a row of pixels is missing, and the text is completely legible. This only impacts SPCs and these are very unlikely to be printed.

6.5 Document changes

| Category | Document | Updated |
|----------------------------|---|---------|
| Deployment | DXS Centrally Hosted SQL Server | |
| | DXS Operating Environment Guidance | |
| | DXS PoC Deployment Detail | |
| | DXS Practice Content Migration | |
| | Technical Survey Overview | |
| | TIO Questionnaire | |
| Knowledge Base | DXS PoC Migration – CSU Instructions | |
| | DXS Program Upgrades and Environment Overview | |
| Technical Overview | PoC Feedback Messaging Process | |
| System Requirements | Release Notes | X |
| | DXS Centrally Hosted SQL Server | |
| | DXS Operating Environment Guidance | |
| | DXS PoC Deployment Detail | |
| | Technical Survey Overview | |
| | TIO Questionnaire | |