



Point-of-Care

Clinical Decision Support

DXS PoC Release Notes Version 5.5.2 – General Release

Date: 3 June 2022

Document Owner: Zayne Glover

Content

1. Product Description	2
2. Releases	2
3. Notes on Deployment	2
4. System Requirements	3
4.1 Recommended Specification	3
4.2 Remote Database Server Specification	3
5. PoC 5.5.2 – General Release	4
5.1 Introduction	5
5.2 Features	5
5.3 Fixes	5
5.4 Known Issues.....	5
5.5 Document changes	6
6. PoC 5.5.1 – Maintenance Release.....	7
6.1 Introduction	8
6.2 Fixes	8
6.3 Known Issues.....	8
6.4 Document changes	8
7. PoC 5.5.1 – General Release	9
7.1 Introduction	10
7.2 Enhancements.....	10
7.3 Fixes	10
7.4 Known Issues.....	10
7.5 Document changes	11

1. Product Description

The DXS Point-of-Care™ product is a content aggregator, which integrates with Electronic Patient Record (EPR) software, ensuring a prime position on the healthcare professional's desktop. This allows content providers to have a direct channel to the healthcare provider during the patient consultations. The DXS product enriches the content with its own relationship-linking process, part of which incorporates the appropriate country-specific diagnostic coding system.

The diagnosis is entered into the EPR and information in DXS Point-of-Care™ relating to that diagnosis is automatically populated for the healthcare provider. The healthcare provider is presented with the specific information related to the patient he or she is consulting, without having to search different databases, which expedites the review.

The DXS Point-of-Care™ product contains clinical data on diseases and procedures as well as patient information on drugs and diseases.

2. Releases

Status	Version	Release Date	Database Schema
Current Release	5.5.2.122	3 June 2022	Version 8 Revision 14
Previous Release	5.5.1.39	13 April 2022	Version 8 Revision 14

3. Notes on Deployment

The upgraded files will be released electronically and will be downloaded with the daily content update. After the update, the upgrade will be pending, and the application will need to restart to deploy the new files.

4. System Requirements

4.1 Recommended Specification

DXS Point-of-Care	Recommended Requirements
Processor	Intel i3 or better
RAM	2GB RAM
Workstation – HDD – Free Space	200MB
Server – HDD – Free Space	4GB (Growth detailed in Operating Environment Guidance document, paragraph 2.5)
Resolution	1280 × 720
Display	32-bit Colour
Operating System	Windows 7, 8, 8.1, 10 supported
Internet Connectivity	Required (HTTP access to uk.dxs-systems.com. HTTP GET and POST Methods allowed)
Internet Speed	5MB (Download), 2MB (Upload) – Please see below for practices that have branch sites

4.2 Remote Database Server Specification

The below Internet speeds are our minimum requirements if a branch site connects remotely via a WAN to the same database as the main site.

DXS Point-of-Care	Recommended Requirements
Processor	Intel i3 or better
RAM	2GB RAM
Workstation – HDD – Free Space	200MB
Server – HDD – Free Space	4GB (Growth detailed in Operating Environment Guidance document, paragraph 2.5)
Resolution	1280 × 720
Display	16 -bit Colour
Operating System	Windows 7, 8, 8.1, 10 supported
Internet Connectivity	Required (HTTP access to uk.dxs-systems.com. HTTP GET and POST Methods allowed)
Internet Speed	10MB (Download), 10MB (Upload)

5. PoC 5.5.2 – General Release

Application	Version	Date	Status
DXS.exe	5.3.2.122	01.06.2022	Updated
DXSSch.exe	5.0.1.22	29.10.2016	–
DXScfg.exe	2.0.0.25	25.02.2005	–
DXSBar.exe	N/A	14.06.2005	–
Launch.exe	2.3.2.122	01.06.2022	Updated
APIProxy.exe	2.5.2.122	01.06.2022	Updated
RegLib.exe	N/A	01.06.2022	Updated
DXS Loader.exe	1.0.0.44	17.03.2021	–
Linker.exe	1.2.2.122	01.06.2022	Updated
HWUpSvc	1.1.0.38	10.02.2021	–

Components	Version	Date	Status
Hwise32x.dll	5.5.2.122	01.06.2022	Updated
Hwise32x.eng	5.5.2.122	01.06.2022	Updated
DXSBar32.dll	N/A		–
DCINPS.dxc	3.3.2.122	01.06.2022	Updated
DCEMISWb.dxc	2.3.2.122	01.06.2022	Updated
DCTPSO.dxc	2.3.2.122	01.06.2022	Updated
Stats.dxp	2.0.2.11	13.08.2013	–
PDPL.dxp	3.0.1.19	11.12.2018	–
Cardiac.dxp	1.0.2.16	10.12.2008	–
ssleay32.dll	1.0.2.15	11.12.2018	–
libeay32.dll	1.0.2.15	11.12.2018	–
FirebirdSql.Data.FirebirdClient.dll	4.6.1.0	06.01.2020	–
HWInput.eng	2.3.2.122	01.06.2022	Updated
LatticeLog.dll	1.4.1.1	06.01.2020	–
LatticeLogBI.dll	1.0.3.1	06.01.2020	–
LatticeLogDAL.dll	1.0.3.1	06.01.2020	–
Microsoft.ServiceBus.dll	3.0.51272.20	06.01.2020	–
Newtonsoft.Json.dll	10.0.1.20720	06.01.2020	–
Lwxspell.dll	1.3.1.0	27.02.2020	–
Editorial.dxp	1.0.2.67	29.07.2021	–

Help Files	Date
DXS.chm	17.05.2007
PDPL.chm	16.07.2009

5.1 Introduction

This version concerns improving user experience. Bug fixes were a crucial focus as well as the complete removal of the Directory of Services module.

5.2 Features

1. Directory of Services module has been removed.
2. Search results that appear under Recent or Favourites are duplicated and no longer replaced.
3. Draft forms are no longer able to save to the patient record.
4. Implementation of a new method of sending mail to Outlook without relying on MAPI.
5. Error telemetry support extended to clinical API calls exceeding 15 seconds.

5.3 Fixes

1. TFS-1201 – Fix for duplicate roles being created for users.
2. TFS-1438 – Fix for the alternative address not displaying on forms.
3. TFS-4164 – Fix for input errors displaying for calculation fields in hidden tables.
4. TFS-4649 – Fix for certain referral links on pathways not working.
5. TFS-4728 – Fix for refreshing of patient data causing PDF view buttons to disappear.
6. TFS-4752 – Fix for form filters being on by default after being turned off.
7. TFS-4800 – Fix for certain Vision sites failing to initialise DXS Point-of-Care.
8. TFS-4865 – Fix for SNOMED codes failing to save for patient education leaflets using TPP.
9. TFS-4902 – Fix for search fields not populating due to mismatch of the coding system.
10. TFS-4905 – Fix for the performance issues experienced at some Vision sites.
11. TFS-4906 – Fix for Favourites not being renamed.
12. TFS-4925 – Fix for exclusive required checkboxes reverting to red colours after saving.
13. TFS-4938 – Fix for excessively large patient records not fully processing and displaying on forms.
14. TFS-4940 – Fix for certain fields intermittently failing to load data.
15. TFS-4963 – Fix for scrolling on forms after clicking on a checkbox.

5.4 Known Issues

1. PRB-106289 – Unable to click on search results.
2. PRB-102104 – Form saved is not the intended form.
3. PRB-112619 – Internal email client is used despite being set to external client.
4. PRB-112622 – Vision sites not updating.

5.5 Document changes

Category	Document	Updated
Deployment	DXS Centrally Hosted SQL Server	
	DXS Operating Environment Guidance	
	DXS PoC Deployment Detail	
	DXS Practice Content Migration	
	Technical Survey Overview	
	TIO Questionnaire	
Knowledge Base	DXS PoC Migration – CSU Instructions	
	DXS Program Upgrades and Environment Overview	
Technical Overview	PoC Feedback Messaging Process	
System Requirements	Release Notes	X
	DXS Centrally Hosted SQL Server	
	DXS Operating Environment Guidance	
	DXS PoC Deployment Detail	
	Technical Survey Overview	
	TIO Questionnaire	

6. PoC 5.5.1 – Maintenance Release

Application	Version	Date	Status
DXS.exe	5.3.1.34	15.03.2022	–
DXSSch.exe	5.0.1.22	29.10.2016	–
DXScfg.exe	2.0.0.25	25.02.2005	–
DXSBar.exe	N/A	14.06.2005	–
Launch.exe	2.3.1.39	05.04.2022	Updated
APIProxy.exe	2.5.1.34	15.03.2022	–
RegLib.exe	N/A	05.04.2022	Updated
DXS Loader.exe	1.0.0.44	17.03.2021	–
Linker.exe	1.2.1.39	05.04.2022	Updated
HWUpSvc	1.1.0.38	10.02.2021	–

Components	Version	Date	Status
Hwise32x.dll	5.5.1.39	05.04.2022	Updated
Hwise32x.eng	5.5.1.39	05.04.2022	Updated
DXSBar32.dll	N/A		–
DCINPS.dxc	3.3.1.34	15.03.2022	–
DCEMISWb.dxc	2.3.1.34	15.03.2022	–
DCTPSO.dxc	2.3.1.34	15.03.2022	–
Stats.dxp	2.0.2.11	13.08.2013	–
PDPL.dxp	3.0.1.19	11.12.2018	–
Cardiac.dxp	1.0.2.16	10.12.2008	–
ssleay32.dll	1.0.2.15	11.12.2018	–
libeay32.dll	1.0.2.15	11.12.2018	–
FirebirdSql.Data.FirebirdClient.dll	4.6.1.0	06.01.2020	–
HWInput.eng	2.3.1.39	05.04.2022	Updated
LatticeLog.dll	1.4.1.1	06.01.2020	–
LatticeLogBI.dll	1.0.3.1	06.01.2020	–
LatticeLogDAL.dll	1.0.3.1	06.01.2020	–
Microsoft.ServiceBus.dll	3.0.51272.20	06.01.2020	–
Newtonsoft.Json.dll	10.0.1.20720	06.01.2020	–
Lwxspell.dll	1.3.1.0	27.02.2020	–
Editorial.dxp	1.0.2.67	29.07.2021	–

Help Files	Date
DXS.chm	17.05.2007
PDPL.chm	16.07.2009

6.1 Introduction

This maintenance release is focussed on a few minor issues that hampered the user's work efficiency.

6.2 Fixes

1. TFS-4793 – Certain emails sent through Outlook fail to send.
2. TFS-4799 – Toolbar search defaults to Education.
3. TFS-4801 – Referral form appears un-editable after changing patients.
4. TFS-4807 – Completed forms section missing one week of data.

6.3 Known Issues

1. PRB-92753 – Duplicate roles listed in the Admin Panel.
2. PRB-106289 – Unable to click on search results.
3. PRB-102104 – Form saved is not the intended form.
4. PRB-112619 – Internal email client is used despite being set to external client.
5. PRB-112621 – Vision sites not initializing due to Isabel.
6. PRB-112622 – Vision sites not updating.

6.4 Document changes

Category	Document	Updated
Deployment	DXS Centrally Hosted SQL Server	
	DXS Operating Environment Guidance	
	DXS PoC Deployment Detail	
	DXS Practice Content Migration	
	Technical Survey Overview	
	TIO Questionnaire	
Knowledge Base	DXS PoC Migration – CSU Instructions	
	DXS Program Upgrades and Environment Overview	
Technical Overview	PoC Feedback Messaging Process	
System Requirements	Release Notes	X
	DXS Centrally Hosted SQL Server	
	DXS Operating Environment Guidance	
	DXS PoC Deployment Detail	
	Technical Survey Overview	
	TIO Questionnaire	

7. PoC 5.5.1 – General Release

Application	Version	Date	Status
DXS.exe	5.3.1.34	15.03.2022	Updated
DXSSch.exe	5.0.1.22	29.10.2016	–
DXScfg.exe	2.0.0.25	25.02.2005	–
DXSBar.exe	N/A	14.06.2005	–
Launch.exe	2.3.1.34	15.03.2022	Updated
APIProxy.exe	2.5.1.34	15.03.2022	Updated
RegLib.exe	N/A	15.03.2022	Updated
DXS Loader.exe	1.0.0.44	17.03.2021	–
Linker.exe	1.2.1.34	15.03.2022	Updated
HWUpSvc	1.1.0.38	10.02.2021	–

Components	Version	Date	Status
Hwise32x.dll	5.5.1.34	15.03.2022	Updated
Hwise32x.eng	5.5.1.34	15.03.2022	Updated
DXSBar32.dll	N/A		–
DCINPS.dxc	3.3.1.34	15.03.2022	Updated
DCEMISWb.dxc	2.3.1.34	15.03.2022	Updated
DCTPSO.dxc	2.3.1.34	15.03.2022	Updated
Stats.dxp	2.0.2.11	13.08.2013	–
PDPL.dxp	3.0.1.19	11.12.2018	–
Cardiac.dxp	1.0.2.16	10.12.2008	–
ssleay32.dll	1.0.2.15	11.12.2018	–
libeay32.dll	1.0.2.15	11.12.2018	–
FirebirdSql.Data.FirebirdClient.dll	4.6.1.0	06.01.2020	–
HWInput.eng	2.3.1.34	15.03.2022	Updated
LatticeLog.dll	1.4.1.1	06.01.2020	–
LatticeLogBI.dll	1.0.3.1	06.01.2020	–
LatticeLogDAL.dll	1.0.3.1	06.01.2020	–
Microsoft.ServiceBus.dll	3.0.51272.20	06.01.2020	–
Newtonsoft.Json.dll	10.0.1.20720	06.01.2020	–
Lwxspell.dll	1.3.1.0	27.02.2020	–
Editorial.dxp	1.0.2.67	29.07.2021	–

Help Files	Date
DXS.chm	17.05.2007
PDPL.chm	16.07.2009

7.1 Introduction

Version 5.5.1 focuses on improving user experience in the 'Forms' section, with various improvements made to aid in filtering and the browsing of completed forms.

7.2 Enhancements

1. A dialog has been added as an additional warning when a patient is changed in the clinical system and a form is loaded with a different patient's information.
2. Information bar is displayed when a filter is applied in the Forms section Inbox folder.
3. Option added to switch between instantly applying form filters and pressing a button to apply the filter.
4. Form Filter options are now limited to valid data available to the user.
5. A dialog message has been added to inform users when no forms are available based on the filters.
6. With the impending change of provider from Patient UK to NHS Choices, users will be asked to amend their Favourites once the former provider's articles are removed.
7. Completed folder has been broken down into weekly folders, including the last 7 days.
8. Non-paginated view has been added back, with the option to switch between it and Paginated view.
9. Additional error tracking has been added to the TPP integration to give more insight into errors received.
10. Directory of Services search results have been removed from the Home screen.

7.3 Fixes

1. TFS-4357 – Fix for form filters causing the product to be extremely slow or crash.
2. TFS-4385 – Fix for completed forms not being archived if they have been deleted.
3. TFS-4387 – Fix for the previous patient's details being retained after changing patients in TPP.
4. TFS-4389 – Fix for search results not opening in the Completed folder.
5. TFS-4634 – Fix for certain pathways not loading because of the CCG branding.
6. TFS-4646 – Fix for column names in the Deleted Forms folder.

7.4 Known Issues

1. PRB-92753 – Duplicate roles listed in the Admin Panel.
2. PRB-106289 – Unable to click on search results.

7.5 Document changes

Category	Document	Updated
Deployment	DXS Centrally Hosted SQL Server	
	DXS Operating Environment Guidance	
	DXS PoC Deployment Detail	
	DXS Practice Content Migration	
	Technical Survey Overview	
	TIO Questionnaire	
Knowledge Base	DXS PoC Migration – CSU Instructions	
	DXS Program Upgrades and Environment Overview	
Technical Overview	PoC Feedback Messaging Process	
System Requirements	Release Notes	X
	DXS Centrally Hosted SQL Server	
	DXS Operating Environment Guidance	
	DXS PoC Deployment Detail	
	Technical Survey Overview	
	TIO Questionnaire	