



Point-of-Care

Clinical Decision Support

DXS PoC Release Notes

Version 5.8.5 – General Release

Date: 9 October 2025

Document Owner: QA Department

Updated by: Zayne Glover

Content

1.	Product Description	2
2.	Releases	2
3.	Notes on Deployment	2
4.	System Requirements	2
5.	PoC 5.8.5 – General Release	3
5.1	Introduction	4
5.2	Features	4
5.3	Fixes	4
5.4	Known Issues.....	5
5.5	Document Changes	5
6.	PoC 5.8.1 – Emergency Release	6
6.1	Introduction	7
6.2	Enhancements.....	7
6.3	Fixes	7
6.4	Known Issues.....	8
6.5	Document Changes	8

1. Product Description

The DXS Point-of-Care™ product is a content aggregator, which integrates with Electronic Patient Record (EPR) software, ensuring a prime position on the healthcare professional's desktop. This allows content providers to have a direct channel to the healthcare provider during the patient consultations. The DXS product enriches the content with its own relationship-linking process, part of which incorporates the appropriate country-specific diagnostic coding system.

The diagnosis is entered into the EPR and information in DXS Point-of-Care™ relating to that diagnosis is automatically populated for the healthcare provider. The healthcare provider is presented with the specific information related to the patient he or she is consulting, without having to search different databases, which expedites the review.

The DXS Point-of-Care™ product contains clinical data on diseases and procedures as well as patient information on drugs and diseases.

2. Releases

Status	Version	Release Date	Database Schema
Current Release	5.8.5.94	6 October 2025	Version 8 Revision 14
Previous Release	5.7.9.27	5 November 2024	Version 8 Revision 14

3. Notes on Deployment

The upgraded files will be released electronically and will be downloaded with the daily content update. After the update, the upgrade will be pending, and the application will need to restart to deploy the new files.

4. System Requirements

Please see “DXS Operating Environment Guidance” document for detailed system specifications:

https://dxsint.sharepoint.com/:w:/g/ETdBSG79zXhAqA1ThgJ0cBoBgae8801SXWI_E1HelBgQrA

5. PoC 5.8.5 – General Release

Application	Version	Date	Status
DXS.exe	5.6.5.94	06.10.2025	Updated
DXSSch.exe	5.1.9.34	25.11.2019	–
DXScfg.exe	2.0.0.25	25.02.2005	–
DXSBar.exe	N/A	14.06.2005	Discontinued
Launch.exe	2.6.5.94	06.10.2025	Updated
APIProxy.exe	2.8.5.94	06.10.2025	Updated
RegLib.exe	N/A	06.10.2025	Updated
DXS Loader.exe	1.0.0.44	17.03.2021	–
Linker.exe	1.5.5.94	06.10.2025	Updated
HWUpSvc.exe	1.5.5.94	06.10.2025	Updated

Components	Version	Date	Status
Hwise32x.dll	5.8.5.94	06.10.2025	Updated
Hwise32x.eng	5.8.5.94	06.10.2025	Updated
DXSBar32.dll	N/A		Discontinued
DCINPS.dxc	3.6.5.94	06.10.2025	Updated
DCEMISWb.dxc	2.6.5.94	06.10.2025	Updated
DCTPSO.dxc	2.6.5.94	06.10.2025	Updated
Stats.dxp	2.0.2.11	13.08.2013	Discontinued
PDPL.dxp	3.0.1.19	11.12.2018	Discontinued
Cardiac.dxp	1.0.2.16	10.12.2008	Discontinued
ssleay32.dll	1.0.2.15	11.12.2018	–
libeay32.dll	1.0.2.15	11.12.2018	–
FirebirdSql.Data.FirebirdClient.dll	4.6.1.0	06.01.2020	–
HWInput.eng	5.8.5.94	06.10.2025	Updated
LatticeLog.dll	1.4.1.1	06.10.2025	Updated
LatticeLogBI.dll	1.0.3.1	06.10.2025	Updated
LatticeLogDAL.dll	1.0.3.1	06.10.2025	Updated
Microsoft.ServiceBus.dll	3.0.51272.20	06.01.2020	–
Newtonsoft.Json.dll	10.0.1.20720	06.01.2020	–
Lwxspell.dll	1.3.1.0	27.02.2020	–
Editorial.dxp	1.0.2.67	29.07.2021	–
ExtAuth.dll	1.2.1.7	06.10.2025	Updated
WebView2Loader.dll	1.0.1054.31	05.02.2024	–

Help Files	Date	Status
DXS.chm	17.05.2007	Discontinued
PDPL.chm	16.07.2009	Discontinued

5.1 Introduction

In version 5.8.5, several key enhancements have been introduced to improve usability and compliance. Form exporting has been updated so that blank forms now include a watermark, and LWD files are exported as PNG rather than PDF files for better accessibility. Additional improvements include displaying “No Data” text as hint text (rather than standard text) in Required and Recommended search fields, new warning prompts when changing Conduit settings or attempting to modify a patient’s NHS number on a referral, and alignment of medicines classifications within Point-of-Care to the latest HDS/BNF standards.

5.2 Features

1. 11045 – Export and print Pathway content as PNG images
2. 10694 – Blank exported forms now include a watermark
3. 5498 – Warning prompt shown when changing Conduit settings
4. 4885 – Required and Recommended search fields display “No Data” as hint text when no results are found
5. 10604 – Support to render Min/Max/Between selection limits for groupers
6. 10773 – Support to render Reusable Custom Input/Error/Complete Validation messages
7. 10922 – Drive space support added for Heartbeat reporting
8. 10074 – Error telemetry logs now include Document ID and SNOMED code for failed patient record attachments
9. 9804 – Medicines classifications in PoC updated to match HDS/BNF standards

5.3 Fixes

1. PRB-138841 – Fixed forms staying in Inbox, while still saving to patient record, instead of moving to Completed
2. 1107 – Fixed hyperlink trigger shifting under hyperlinked text in vertically centred cells
3. 6311 – Multiple users editing a form can result in changes being overwritten
4. 7882 – Resolved issue with several SNOMED passback codes not saving in TPP
5. 10989 – Hidden forms no longer show in the Summary pane when intended to hide
6. 10536 – Validation message now remains visible when cancelling “Complete” to edit recommended fields
7. 10867 – Fixed form attachment failure caused by an extra “N” added to an NHS number
8. 10873 – Removed incorrect “Major Problems” text from the Severe Problems tag
9. 10926 – PDF metadata is now correctly embedded when a form is emailed
10. 10977 – Signature field now clears correctly when clearing a form
11. 10992 – Removed e-RS submission text from document actions while e-RS is unsupported
12. 11100 – Fixed PathBlock alignment issues in Suffolk ICB content
13. 11124 – Resolved crash when loading forms containing infinite recursion
14. 11266 – Double-clicking search button/enter no longer stalls pending and future searches
15. 11367 – Field input tips no longer appear through system dialogues or external apps
16. 11402 – Editing of the NHS number in a populated form is now prevented
17. 10839 – Removed outdated “?” button from the search bar
18. 11467 – Fixed Product listing being removed from Add/Remove programs

5.4 Known Issues

1. PRB-86507 – Duplicate sites appearing in the database
2. PRB-92753 – Duplicate Roles listed in Admin panel
3. PRB-121017 – Print spooler causes product to hang when saving to the patient record
4. PRB-126423 – Bell notification for assigned forms not appearing over VPN connections
5. PRB-141806 – File certificates missing after upgrades
6. PRB-147345 – TPP API multiple app handling issue
7. 6999 – Save to Disk option changes white text to black text
8. 7541 – Placeholder field colour properties do not display correctly
9. 8545 – Form titles can be cut off due to a character limit
10. 9963 – Some summary hints are cut off when displayed
11. 10770 – Strict Age Validation can be bypassed
12. 10899 – PathEdit PDF export alignment issues
13. 10917 – Ordering of clinical entries added on the same day can be inconsistent
14. 10919 – TPP clinical entry without a SNOMED code displays a broken or missing value
15. 11241 – Second Windows user inherits first user's login
16. 11453 – Custom cell colour appears in auto-populated cells on PDF export

5.5 Document Changes

Category	Document	Updated
Deployment	DXS Centrally Hosted SQL Server	
	DXS PoC Deployment Detail	
	DXS Practice Content Migration	
	Technical Survey Overview	
	TIO Questionnaire	
Knowledge Base	DXS PoC Migration – CSU Instructions	
	DXS Program Upgrades and Environment Overview	
Technical Overview	PoC Feedback Messaging Process	
System Requirements	Release Notes	X
	DXS Operating Environment Guidance	

6. PoC 5.8.1 – Emergency Release

Application	Version	Date	Status
DXS.exe	5.6.1.55	31.03.2025	Updated
DXSSch.exe	5.1.9.34	25.11.2019	–
DXScfg.exe	2.0.0.25	25.02.2005	–
DXSBar.exe	N/A	14.06.2005	Discontinued
Launch.exe	2.6.1.55	31.03.2025	Updated
APIProxy.exe	2.8.1.55	31.03.2025	Updated
RegLib.exe	N/A	31.03.2025	Updated
DXS Loader.exe	1.0.0.44	17.03.2021	–
Linker.exe	1.5.1.55	31.03.2025	Updated
HWUpSvc.exe	1.5.1.55	31.03.2025	Updated

Components	Version	Date	Status
Hwise32x.dll	5.8.1.55	31.03.2025	Updated
Hwise32x.eng	5.8.1.55	31.03.2025	Updated
DXSBar32.dll	N/A		Discontinued
DCINPS.dxc	3.6.1.55	31.03.2025	Updated
DCEMISWb.dxc	2.6.1.55	31.03.2025	Updated
DCTPSO.dxc	2.6.1.55	31.03.2025	Updated
Stats.dxp	2.0.2.11	13.08.2013	Discontinued
PDPL.dxp	3.0.1.19	11.12.2018	Discontinued
Cardiac.dxp	1.0.2.16	10.12.2008	Discontinued
ssleay32.dll	1.0.2.15	11.12.2018	–
libeay32.dll	1.0.2.15	11.12.2018	–
FirebirdSql.Data.FirebirdClient.dll	4.6.1.0	06.01.2020	–
HWInput.eng	5.8.1.55	31.03.2025	Updated
LatticeLog.dll	1.4.1.1	06.01.2020	–
LatticeLogBI.dll	1.0.3.1	06.01.2020	–
LatticeLogDAL.dll	1.0.3.1	06.01.2020	–
Microsoft.ServiceBus.dll	3.0.51272.20	06.01.2020	–
Newtonsoft.Json.dll	10.0.1.20720	06.01.2020	–
Lwxspell.dll	1.3.1.0	27.02.2020	–
Editorial.dxp	1.0.2.67	29.07.2021	–
ExtAuth.dll	1.2.1.7	31.03.2025	Updated
WebView2Loader.dll	1.0.1054.31	05.02.2024	–

Help Files	Date	Status
DXS.chm	17.05.2007	Discontinued
PDPL.chm	16.07.2009	Discontinued

6.1 Introduction

Version 5.8.1 introduces a new “Hidden” referral form behaviour, allowing certain forms to be accessible only via a direct link. This ensures more controlled access where necessary. Additionally, the form metadata versioning has been updated to prevent older product versions from displaying forms that include unsupported new features, helping maintain compatibility and reduce errors. As development on e-RS continues, it has been temporarily disabled in Version 5.8.1.

6.2 Enhancements

1. e-RS has been disabled for further development.
2. Added error telemetry for failed Signature checks during start-up.
3. An override has been implemented to address failures during the Signature check at start-up.
4. The metadata version of documents has been updated to prevent incompatible functionality from not working on older versions of PoC.
5. Added “Hidden” referral form behaviour where certain forms are only accessible via a direct link from a relevant criteria page.
6. The ability to tag documents with specific metadata that is included in PDF exports has been added.

6.3 Fixes

1. PRB-106289 – Fix for Search results freezing.
2. 9915 – Fix for checkbox Input Error icon displaying incorrectly.
3. 10025 – Fix for numeric fields accepting invalid values and units duplicating.
4. 10036 – Fix for switching patients and EMIS instances resulting in a form bypassing validation.
5. 10164 – Fix for invalid date values changing when navigating away and back to the same form.
6. 10196 – Fix for the date disambiguation pop-up no longer occurring.
7. 10207 – Fix for forms summary actions not exporting Reasonable Adjustments notification banner.
8. 10209 – Fix for the error codes missing on error messages.
9. 10358 – Fix for default input on Radio, Check, and Tick Lookup lists not reflecting on PDF export.
10. 10472 – Fix for PDF exports failing to render via document actions on some forms.

6.4 Known Issues

1. PRB-86507 – Duplicate sites appearing in the database.
2. PRB-92753 – Duplicate Roles listed in Admin panel.
3. PRB-121017 – Print spooler causes product to hang on saving to the patient record.
4. PRB-126423 – Bell notification for assigned forms not appearing over VPN connections.
5. PRB-141806 – File certificates missing after upgrades.
6. 1107 – Hyperlink trigger shifts under hyperlinked text in vertically centered cells.
7. 5376 – Warning for locked forms only showing after editing when attempting to save.
8. 7541 – Placeholder field colour properties do not display correctly.
9. 9963 – Some summary hints are cut off when displayed.
10. 10089 – The e-RS Only filter option in the Forms section disappears intermittently.
11. 10460 – e-RS cover letter not saving back for test patient that has special characters within their name and surname.

6.5 Document Changes

Category	Document	Updated
Deployment	DXS Centrally Hosted SQL Server	
	DXS PoC Deployment Detail	
	DXS Practice Content Migration	
	Technical Survey Overview	
	TIO Questionnaire	
Knowledge Base	DXS PoC Migration – CSU Instructions	
	DXS Program Upgrades and Environment Overview	
Technical Overview	PoC Feedback Messaging Process	
System Requirements	Release Notes	X
	DXS Operating Environment Guidance	