



# Point-of-Care

Clinical Decision Support

## DXS PoC Release Notes

Customer-Facing

Version 5.9.0 – General Release

Date: 23 March 2026

Document Owner: QA Department

Updated by: Zayne Glover

# Content

1. Product Description .....	2
2. Releases .....	2
3. Notes on Deployment .....	2
4. System Requirements .....	2
5. PoC 5.9.0 – General Release .....	3
5.1 Introduction .....	3
5.2 Features .....	3
5.3 Fixes .....	3
5.4 Known Issues.....	3
5.5 Known Third-Party Limitations .....	4

## 1. Product Description

DXS Point-of-Care is a clinical decision support tool designed to help healthcare professionals make faster, safer, and more accurate decisions at the point of care. It provides access to automatically updated clinical pathways, referral forms, patient information leaflets, and medicines information – all in one place.

## 2. Releases

Status	Version	Release Date	Database Schema
Current Release	5.9.0.61	20 March 2026	Version 8 Revision 14
Previous Release	5.8.5.106	17 December 2025	Version 8 Revision 14

## 3. Notes on Deployment

The upgraded files will be downloaded automatically as part of the daily content update.

Once the update is complete, the upgrade will appear as *pending*. To deploy the new files, the application must be restarted. After restarting, DXS Point-of-Care will be updated to the latest version.

If you experience any issues with the automated upgrade process, please contact [Support](#).

## 4. System Requirements

For detailed system specifications, please refer to the *DXS Operating Environment Guidance* document:

<https://storage.googleapis.com/images003/DXSDocuments/DXS%20Operating%20Environment%20Guidance.pdf>

## 5. PoC 5.9.0 – General Release

### 5.1 Introduction

Version 5.9.0.61 introduces enhancements to email authentication, system security, and overall usability. The primary focus of this release is the implementation of OAuth 2.0 email authentication, enabling modern, token-based authentication aligned with current security standards.

This release also strengthens secure communications, improves configuration flexibility, and enhances system reliability through updated telemetry and retry mechanisms.

### 5.2 Features

1. **Azure-10820** – Added the ability to exclude selected fields during form completion.
2. **Azure-11081** – Introduced support for OAuth 2.0 email authentication.
3. **Azure-11936** – Updated error telemetry to use Transport Layer Security (TLS) 1.2 for enhanced security.
4. **Azure-11937** – Implemented automated updates for File Signature Thumbprints.
5. **Azure-11947** – Added retry capability for File Signature Thumbprint verification checks.

### 5.3 Fixes

1. **Azure-1093** – Fixed an issue where the cursor incorrectly hovered over the Search hint text in the toolbar.
2. **Azure-10878** – Resolved an issue where forms saved to the patient record remained in the PoC Inbox instead of moving to the Completed folder.
3. **Azure-10899** – Fixed alignment issues when exporting Pathways to PDF.
4. **Azure-10919** – Fixed an issue where data from TPP displayed incorrectly when no SNOMED code was present.
5. **Azure-11753** – Resolved an issue where special characters caused PDF metadata to fail during export.
6. **Azure-11808** – Restored the ability to export form summaries from dated Completed Form folders.
7. **Azure-11877** – Fixed an issue where validation descriptions were not used as fallback messages for form validation errors.
8. **Azure-12069** – Fixed an issue where shared group rule error messages displayed the incorrect values.
9. **Azure-12076** – Resolved an issue where EMIS filtering could cause the toolbar to become unresponsive.
10. **Azure-12142** – Fixed an issue where emails sent via the default DXS email configuration could bounce back.

### 5.4 Known Issues

1. **PRB-92753** – Duplicate roles listed in the Admin panel.  
**Workaround:** Please contact Support to remove duplicate roles.

2. **PRB-121017** – Application may hang when saving to the patient record due to print spooler behaviour.  
**Workaround:**  
Please contact Support, or apply the following configuration change:  
Navigate to:  
HKEY\_LOCAL\_MACHINE\SYSTEM\CurrentControlSet\Control\Terminal Server\Sysprocs
  - Create a new **DWORD** value named `splwow64.exe`.
  - Set the value data to **0** to enable automatic termination.
3. **PRB-147345** – TPP API timeout when loading patient demographics onto forms.  
**Workaround:** This issue is currently under investigation. Please contact Support.
4. **Azure-10917** – Clinical entries added on the same day may appear in an inconsistent order due to lack of timestamping.
5. **Azure-12022** – Admin draft override may allow actions on forms that are already completed or deleted.
6. **Azure-12024** – Forms may be deleted from the Summary section while locked and still in use.
7. **Azure-12050** – Favourited documents in the Forms section display the name used when favourited, rather than the original form name, which may cause confusion when assigned to another user.
8. **Azure-12052** – Validation errors may not display on form load and only appear after interacting with a field.

## 5.5 Known Third-Party Limitations

1. **PRB-126423** – Bell notifications for assigned forms may not appear when using VPN connections.
2. **Azure-11241** – In TPP integrations, a second Windows user may inherit the first user's login session.
3. **Azure-11996** – Centile codes in EMIS Growth Charts may prevent Height and Weight codes from populating in search fields.